


# LMK EyeFit

LUMINANCE & COLOR  
CHARACTERIZATION  
FOR AR/VR DEVICES



**LMK6-30**

**LMK6-12**

| NED lens type           | TT NED 40                               | TT NED 40 Autofocus                     |
|-------------------------|---|---|
| Focus operation         | Manually controlled focus               | Electronically controlled focus         |
| Focal power             | 0 dpt. to +2 dpt.                       | -5 dpt. to +10 dpt.                     |
| Entrance pupil position | First optical surface                   | First optical surface                   |
| Entrance pupil diameter | 2-6 mm<br>(<2.5 mm diffraction limited) | 2-6 mm<br>(<2.5 mm diffraction limited) |
| FOV (total)             | 30° (H) × 26.4° (V)                     | 30° (H) × 26.2° (V)                     |
| FOV (total)             | 21.2° (H) × 14.8° (V)                   | 21.2° (H) × 14.6° (V)                   |